





**Bio One** Binocular Microscope is a great value for an introductory microscope. The **Bio One** Series comes standard with 10x/18mm eyepieces, **Achromat** 4,10,40,100x objectives, mechanical stage with specimen holder, abbe 1.25 condenser with iris, 6v/20w illumination, CE approved electronics, dust cover, blue filter, immersion oil, 5 year mechanical warranty and 1 year electrical warranty

**Bio View** Binocular or Trinocular Microscope offers great optics and mechanics at an affordable price. The **Bio View** Series comes standard with 10x/18mm eyepieces, **Infinity Corrected Achromat** 4,10,40,100x objectives, 142mmx135mm mechanical stage with specimen holder, abbe 1.25 condenser with iris(marked for objective settings), 6v/20w illumination, CE approved electronics, dust cover, blue filter, immersion oil, 5 year mechanical warranty and 1 year electrical warranty

**Bio Scope** Binocular or Trinocular Microscope offers superb optics and mechanics at a competitive price. The **Bio Scope** Series comes standard with 10x/20mm eyepieces, **Infinity Corrected Plan** 4,10,40,100x objectives, 156mmx138mm mechanical stage with specimen holder, abbe 1.25 condenser with iris(marked for objective settings), Koehler illumination with a 6v/30w halogen bulb, CE approved electronics, dust cover, blue filter, immersion oil, 5 year mechanical warranty and 1 year electrical warranty



**Bio Scope II** Binocular or Trinocular Microscope with beautiful optics and mechanics at a great price. The **Bio Scope II** Series comes standard with 10x/22mm eyepieces, **Infinity Corrected Plan** 4,10,40,100x objectives, Anodized 216mmx150mm mechanical stage with specimen holder, abbe 1.25 condenser with iris (marked for objective settings), Koehler illumination with a 3w LED bulb, CE approved electronics, dust cover, blue filter, immersion oil, 5 year mechanical warranty and 1 year electrical warranty



Shown with Tilting Head and Ceramic Stage



Shown with Trinocular Head

**Bio Imager** Microscope offers outstanding optics and mechanics at an affordable price. The **Bio Imager** Series comes standard with 10x/22mm eyepieces, **Infinity Corrected Plan** 4,10,40,100x objectives, 218mmx150mm mechanical stage with specimen holder, abbe 1.25 condenser with iris, Koehler illumination with 6v/30w bulb, **Contrast Enhancement Filter** (not your standard blue filter), CE approved electronics, dust cover, immersion oil, 5 year warranty

**Optional Items:** Ceramic Coated Mechanical Stage, Tilting Ergonomic Head, Trinocular Bridge, Multi-Observation, 2x, 20x, 50x Oil, 60x Infinity Plan Objectives, Fluorite Optics, Fluorescents, Polarizer/Analyzer



Shown with Binocular Head

< **Bio Lab** Trinocular Microscope offers excellent optics and mechanics at an economical price. The **Bio Lab** Series comes standard with 10x/20mm eyepieces, **Infinity Corrected Plan** 4,10,40,100x objectives, 188mmx150mm ceramic coated mechanical stage with specimen holder, abbe 1.25 condenser with iris(marked for objective settings), 12v/20w illumination, CE approved electronics, dust cover, blue filter, immersion oil, 5 year mechanical warranty and 1 year electrical warranty



**Bio Imager II** Microscope offers bet-in-class performance and value for clinical and research applications. Newly designed NIS infinity Plan optics provide crisp sharp images. The **Bio Imager II** comes standard with 10x/22mm eyepieces, **Infinity Corrected Plan** 4,10,40,100x objectives, 302mmx152mm mechanical stage with Gorilla Glass insert, specimen holder and X/Y friction adjustable knobs, abbe 1.25 condenser with iris, Koehler illumination with LED 25,000 hrs. bulb, CE approved electronics, dust cover, immersion oil, 5 year warranty, 1 year LED bulb warranty

**Optional Items:** Tilting Ergonomic Head, Trinocular Bridge, Multi-Observation, 2x, 20x, 50x Oil, 60x Infinity Plan Objectives, APO Optics, Fluorescents, Polarizer/Analyzer

